

IMPLANT CATALOGUE

EFFECTIVE JANUARY 2015




MENTOR[®]
BREAST IMPLANTS



APPROVED FOR SWITZERLAND

© Mentor Worldwide LLC 2014 025632-141126-CH-001/15

In 2002, we completed a new plant in Leiden, The Netherlands. At 6,000 square meters, it is the largest of its kind in the world. This plant contains an automated breast implant dipping machine, which produces every implant shell in a precisely uniform manner. To meet an ever-increasing demand for our products, the new plant produces more than 300,000 breast implants a year, distributed to more than 100 nations. At Mentor all our products are made under strict standards of design and testing.





Founded in 1969, Mentor Worldwide LLC is a leading manufacturer and supplier of medical products for the global Aesthetic medicine market and Reconstructive surgery. As the world's leading maker of high-quality breast implants for over 30 years, our experience results in quality products that you can rely on.

Headquartered in Santa Barbara, California, Mentor has manufacturing and research operations in the United States and the Netherlands. In 1995, Mentor became the first breast implant manufacturer to achieve CE-marking for its products.

Mentor has worked closely with surgeons throughout the world to develop innovative products such as CONTOUR PROFILE™ Tissue Expanders, BECKER™ Expander/Breast Implants, and the CPG™ Gel Breast Implants.

Mentor's manufacturing capability is exceptional. We have substantial experience in biomaterials and biopolymers manufacturing. At the same time, our regulatory expertise has enabled the Company to bring products to market in the world's most stringently regulated health care environments.

Mentor is a manufacturer with a strong customer focus. Our sales representatives are skilled in technical product issues, enabling them to act as a valuable information resource. For questions or additional information, please contact your local Mentor branch or Johnson & Johnson office (listed on back cover) or distributor.

On Jan 23, 2009, Mentor was acquired by Johnson & Johnson and is a stand-alone business unit reporting through ETHICON Inc., a leading provider of suture, mesh, hemostats and other products for a wide range of surgical procedures. By joining forces, we aspire to be the trusted global leader in aesthetic medicine among both consumers and professionals by providing a broad range of innovative, science and clinical-based solutions to maintain, enhance, and restore self-esteem and quality of life. Simply said, our goal is to make life more beautiful.



CONTENTS

MemoryGel™ Breast Implants,
A More Evolved Implant 1

About MemoryGel™ Implants
MemoryGel™ Implants Product Scale

All Cohesive, All the Time™ 2

Round Gel Breast Implants, Cohesive I™ 3

Smooth Round Breast Implants, Moderate Profile, Cohesive I™
Smooth Round Breast Implants, Moderate Plus Profile, Cohesive I™
Smooth Round Breast Implants, High Profile, Cohesive I™
Smooth Round Breast Implants, Ultra High Profile, Cohesive I™
SILTEX™ Round Breast Implants, Moderate Profile, Cohesive I™
SILTEX™ Round Breast Implants, Moderate Plus Profile, Cohesive I™
SILTEX™ Round Breast Implants, High Profile, Cohesive I™
SILTEX™ Round Breast Implants, Ultra High Profile, Cohesive I™

Round Gel Breast Implants, Cohesive II™ 7

SILTEX™ Round Breast Implants, Moderate Profile, Cohesive II™
SILTEX™ Round Breast Implants, Moderate Plus Profile, Cohesive II™
SILTEX™ Round Breast Implants, High Profile, Cohesive II™

CPG™ Gel Breast Implants, Cohesive III™ 9

CPG™ Implants 331, Cohesive III™, Tall Height, Moderate Projection
CPG™ Implants 321, Cohesive III™, Medium Height, Moderate Projection
CPG™ Implants 311, Cohesive III™, Low Height, Moderate Projection
CPG™ Implants 332, Cohesive III™, Tall Height, Moderate Plus Projection
CPG™ Implants 322, Cohesive III™, Medium Height, Moderate Plus Projection
CPG™ Implants 312, Cohesive III™, Low Height, Moderate Plus Projection
CPG™ Implants 333, Cohesive III™, Tall Height, High Projection
CPG™ Implants 323, Cohesive III™, Medium Height, High Projection
CPG™ Implants 313, Cohesive III™, Low Height, High Projection
CONTOUR PROFILE™ Tissue Expanders / CPG™ Gel Breast Implants
Sizing Guide – Two-Stage Breast Reconstruction

BECKER™ Expander/Breast Implants 15

SILTEX™ Round BECKER™ 25 Expander/Breast Implants, Cohesive I™
SILTEX™ Round BECKER™ 50 Expander/Breast Implants, Cohesive I™
SILTEX™ CONTOUR PROFILE™ BECKER™ 35 Expander/Breast Implants,
Cohesive II™

Sterile Resterilizable Gel Breast Implant Sizers 16

Round Breast Implant Sizers, Moderate Profile
Round Breast Implant Sizers, Moderate Plus Profile
Round Breast Implant Sizers, High Profile
Round Breast Implant Sizers, Ultra High Profile

CPG™ Gel Breast Implant Sizers.....	18
CPG™ Implant Sizers 331, Cohesive III™, Tall Height, Moderate Projection	
CPG™ Implant Sizers 321, Cohesive III™, Medium Height, Moderate Projection	
CPG™ Implant Sizers 311, Cohesive III™, Low Height, Moderate Projection	
CPG™ Implant Sizers 332, Cohesive III™, Tall Height, Moderate Plus Projection	
CPG™ Implant Sizers 322, Cohesive III™, Medium Height, Moderate Plus Projection	
CPG™ Implant Sizers 312, Cohesive III™, Low Height, Moderate Plus Projection	
CPG™ Implant Sizers 333, Cohesive III™, Tall Height, High Projection	
CPG™ Implant Sizers 323, Cohesive III™, Medium Height, High Projection	
CPG™ Implant Sizers 313, Cohesive III™, Low Height, High Projection	
CPX™ Breast Tissue Expanders	21
CONTOUR PROFILE™ Expander, Low Height, Styles 8100, 9100	
CONTOUR PROFILE™ Expander, Medium Height, Styles 8200, 9200	
CONTOUR PROFILE™ Expander, Tall Height, Styles 8300, 9300	
Other Tissue Expanders.....	23
Smooth Round Tissue Expanders	
Smooth Elliptical Tissue Expanders	
Smooth Crescent Tissue Expanders	
Smooth Rectangle Tissue Expanders	
Sterile Saline Breast Implant Sizers	24
Round Breast Implant Sizers, Moderate Profile	
Round Breast Implant Sizers, Moderate Plus Profile	
Round Breast Implant Sizers, High Profile	
CONTOUR PROFILE™ Saline Breast Implant Sizers, High Profile	
Accessories	26
Injection Port Detector	
Injection Dome	
Important Safety Information.....	27
Notes.....	28
Contact Information	29

MemoryGel™

BREAST IMPLANTS

A More Evolved Implant

- Strict US and EU manufacturing standards
- Highly compliant shell for ease of placement
- Lifetime product replacement policy
- PatientSafe coverage limited warranty or PatientCare Promise limited warranty†
- Moderate, Moderate Plus, High, and Ultra High projections
- Proprietary silicone gel formula for a natural touch that resembles breast tissue

About MemoryGel™ Breast Implants

MENTOR® MemoryGel™ Breast Implants feature a proprietary cohesive silicone gel formulation used to fill all MENTOR® Silicone Gel-Filled Breast Implants around the world. By varying the components of the gel, Mentor is able to produce a wide selection of products ranging from very soft to very firm.

MemoryGel™ Breast Implants Product Scale

Product Profile	Cohesive I™	Cohesive II™	Cohesive III™
Round Moderate Plus	X	X	
Round Moderate	X	X	
Round High	X	X	
Round Ultra High	X		
CPG™ Gel Breast Implant			X
Round BECKER™ 25 Expander/Breast Implant	X		
Round BECKER™ 50 Expander/Breast Implant	X		
CONTOUR PROFILE™ BECKER™ 35 Expander/Breast Implant		X	

† In case of eligibility, depending on the date of surgery, the implant will fall under one of the above warranties.

All Cohesive, All the Time.™



Advanced MemoryGel™ Breast Implants.
The more evolved implant.

All MENTOR® Silicone Gel-Filled Breast Implants contain gel that is cohesive, safe and aesthetically pleasing.

Cohesive I™

The standard cohesive level gel used in MENTOR® Breast Implants. This is the softest gel and has been preferred for years in textured and smooth Round Moderate, Moderate Plus, High, and Ultra High Profile Gel Breast Implants as well as the BECKER™ Expander/Breast Implants.

Cohesive II™

A slightly firmer gel, for surgeons wanting a firmer feeling implant. This gel is used in textured Round Moderate, Moderate Plus, High Profile Gel Breast Implants and CONTOUR PROFILE™ BECKER™ Expander/Breast Implants.

Cohesive III™

Mentor's most cohesive gel, providing shape retention with a pleasing level of firmness for optimal aesthetic results. This gel is used in the MENTOR® CPG™ Gel Breast Implants.

MEMORYGEL™ BREAST IMPLANTS

Round Gel Breast Implants, Cohesive I™ CONT'D

Smooth Round Breast Implants, Moderate Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	10,0 cm	2,2 cm	350-7125 BC	RSZ-7125
150 cc	10,3 cm	2,3 cm	350-7150 BC	RSZ-7150
175 cc	11,2 cm	2,4 cm	350-7175 BC	RSZ-7175
200 cc	11,7 cm	2,5 cm	350-7200 BC	RSZ-7200
225 cc	12,2 cm	2,6 cm	350-7225 BC	RSZ-7225
250 cc	12,3 cm	2,8 cm	350-7250 BC	RSZ-7250
275 cc	13,2 cm	2,9 cm	350-7275 BC	RSZ-7275
300 cc	13,5 cm	3,0 cm	350-7300 BC	RSZ-7300
325 cc	13,9 cm	3,0 cm	350-7325 BC	RSZ-7325
350 cc	14,2 cm	3,1 cm	350-7350 BC	RSZ-7350
375 cc	14,4 cm	3,2 cm	350-7375 BC	RSZ-7375
400 cc	14,5 cm	3,2 cm	350-7400 BC	RSZ-7400
450 cc	14,9 cm	3,2 cm	350-7450 BC	RSZ-7450
500 cc	15,2 cm	3,4 cm	350-7500 BC	RSZ-7500
550 cc	15,9 cm	3,6 cm	350-7550 BC	RSZ-7550
600 cc	16,5 cm	3,7 cm	350-7600 BC	RSZ-7600
700 cc	17,4 cm	3,9 cm	350-7700 BC	RSZ-7700
800 cc	18,2 cm	4,1 cm	350-7800 BC	RSZ-7800

Smooth Round Breast Implants, Moderate Plus Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	8,9 cm	2,8 cm	350-1251 BC	RSZ-1251
150 cc	9,5 cm	2,9 cm	350-1501 BC	RSZ-1501
175 cc	10,0 cm	3,1 cm	350-1751 BC	RSZ-1751
200 cc	10,5 cm	3,2 cm	350-2001 BC	RSZ-2001
225 cc	10,9 cm	3,3 cm	350-2251 BC	RSZ-2251
250 cc	11,3 cm	3,4 cm	350-2501 BC	RSZ-2501
275 cc	11,7 cm	3,5 cm	350-2751 BC	RSZ-2751
300 cc	12,0 cm	3,6 cm	350-3001 BC	RSZ-3001
325 cc	12,3 cm	3,8 cm	350-3251 BC	RSZ-3251
350 cc	12,5 cm	3,9 cm	350-3501 BC	RSZ-3501
375 cc	12,8 cm	4,0 cm	350-3751 BC	RSZ-3751
400 cc	13,1 cm	4,0 cm	350-4001 BC	RSZ-4001
450 cc	13,6 cm	4,2 cm	350-4501 BC	RSZ-4501
500 cc	14,1 cm	4,3 cm	350-5001 BC	RSZ-5001
550 cc	14,6 cm	4,5 cm	350-5501 BC	RSZ-5501
600 cc	15,0 cm	4,6 cm	350-6001 BC	RSZ-6001
700 cc	15,8 cm	4,9 cm	350-7001 BC	RSZ-7001
800 cc	16,5 cm	5,1 cm	350-8001 BC	RSZ-8001

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

Round Gel Breast Implants, Cohesive I™ CONT'D

Smooth Round Breast Implants, High Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	8,3 cm	3,5 cm	350-1254 BC	RSZ-1254
150 cc	8,8 cm	3,7 cm	350-1504 BC	RSZ-1504
175 cc	9,3 cm	3,9 cm	350-1754 BC	RSZ-1754
200 cc	9,7 cm	4,0 cm	350-2004 BC	RSZ-2004
225 cc	10,1 cm	4,2 cm	350-2254 BC	RSZ-2254
250 cc	10,5 cm	4,3 cm	350-2504 BC	RSZ-2504
275 cc	10,8 cm	4,4 cm	350-2754 BC	RSZ-2754
300 cc	11,1 cm	4,5 cm	350-3004 BC	RSZ-3004
325 cc	11,4 cm	4,6 cm	350-3254 BC	RSZ-3254
350 cc	11,7 cm	4,8 cm	350-3504 BC	RSZ-3504
375 cc	12,0 cm	4,8 cm	350-3754 BC	RSZ-3754
400 cc	12,2 cm	5,0 cm	350-4004 BC	RSZ-4004
425 cc	12,5 cm	5,0 cm	350-4254 BC	RSZ-4254
450 cc	12,8 cm	5,1 cm	350-4504 BC	RSZ-4504
500 cc	13,2 cm	5,3 cm	350-5004 BC	RSZ-5004
550 cc	13,6 cm	5,5 cm	350-5504 BC	RSZ-5504
600 cc	14,0 cm	5,6 cm	350-6004 BC	RSZ-6004
650 cc	14,4 cm	5,7 cm	350-6504 BC	RSZ-6504
700 cc	14,8 cm	5,8 cm	350-7004 BC	RSZ-7004
800 cc	15,5 cm	6,0 cm	350-8004 BC	RSZ-8004

Smooth Round Breast Implants, Ultra High Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
135 cc	7,8 cm	4,1 cm	350-5135 BC	RSZ-5135
160 cc	8,2 cm	4,3 cm	350-5160 BC	RSZ-5160
185 cc	8,4 cm	4,4 cm	350-5185 BC	RSZ-5185
215 cc	8,7 cm	4,5 cm	350-5215 BC	RSZ-5215
240 cc	9,0 cm	4,7 cm	350-5240 BC	RSZ-5240
270 cc	9,3 cm	4,8 cm	350-5270 BC	RSZ-5270
295 cc	9,7 cm	5,0 cm	350-5295 BC	RSZ-5295
320 cc	9,8 cm	5,1 cm	350-5320 BC	RSZ-5320
350 cc	10,1 cm	5,2 cm	350-5350 BC	RSZ-5350
375 cc	10,4 cm	5,3 cm	350-5375 BC	RSZ-5375
400 cc	10,6 cm	5,4 cm	350-5400 BC	RSZ-5400
430 cc	10,9 cm	5,6 cm	350-5430 BC	RSZ-5430
455 cc	11,2 cm	5,7 cm	350-5455 BC	RSZ-5455
480 cc	11,4 cm	5,8 cm	350-5480 BC	RSZ-5480
535 cc	12,0 cm	6,1 cm	350-5535 BC	RSZ-5535
590 cc	12,5 cm	6,3 cm	350-5590 BC	RSZ-5590

*Resterilizable up to 10 times.

MEMORYGEL™ BREAST IMPLANTS

Round Gel Breast Implants, Cohesive I™

SILTEX™ Round Breast Implants, Moderate Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
100 cc	8,8 cm	2,5 cm	354-1007	RSZ-7100
125 cc	9,3 cm	2,8 cm	354-1257	RSZ-7125
150 cc	10,2 cm	2,7 cm	354-1507	RSZ-7150
175 cc	10,7 cm	2,8 cm	354-1757	RSZ-7175
200 cc	11,2 cm	2,8 cm	354-2007	RSZ-7200
225 cc	11,4 cm	3,0 cm	354-2257	RSZ-7225
250 cc	11,5 cm	3,3 cm	354-2507	RSZ-7250
275 cc	12,4 cm	3,4 cm	354-2757	RSZ-7275
300 cc	12,6 cm	3,5 cm	354-3007	RSZ-7300
325 cc	12,9 cm	3,6 cm	354-3257	RSZ-7325
350 cc	13,4 cm	3,7 cm	354-3507	RSZ-7350
375 cc	13,4 cm	3,8 cm	354-3757	RSZ-7375
400 cc	13,5 cm	3,9 cm	354-4007	RSZ-7400
450 cc	13,9 cm	4,1 cm	354-4507	RSZ-7450
500 cc	14,2 cm	4,3 cm	354-5007	RSZ-7500
550 cc	14,8 cm	4,4 cm	354-5507	RSZ-7550
600 cc	15,4 cm	4,5 cm	354-6007	RSZ-7600
700 cc	16,8 cm	4,3 cm	354-7007	RSZ-7700
800 cc	17,2 cm	4,6 cm	354-8007	RSZ-7800

SILTEX™ Round Breast Implants, Moderate Plus Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
100 cc	8,1 cm	2,7 cm	354-1001	RSZ-1001
125 cc	8,8 cm	2,9 cm	354-1251	RSZ-1251
150 cc	9,4 cm	3,0 cm	354-1501	RSZ-1501
175 cc	10,0 cm	3,2 cm	354-1751	RSZ-1751
200 cc	10,5 cm	3,3 cm	354-2001	RSZ-2001
225 cc	10,9 cm	3,5 cm	354-2251	RSZ-2251
250 cc	11,3 cm	3,6 cm	354-2501	RSZ-2501
275 cc	11,7 cm	3,7 cm	354-2751	RSZ-2751
300 cc	12,0 cm	3,7 cm	354-3001	RSZ-3001
325 cc	12,3 cm	3,8 cm	354-3251	RSZ-3251
350 cc	12,6 cm	3,8 cm	354-3501	RSZ-3501
375 cc	12,9 cm	3,9 cm	354-3751	RSZ-3751
400 cc	13,2 cm	4,0 cm	354-4001	RSZ-4001
450 cc	13,7 cm	4,1 cm	354-4501	RSZ-4501
500 cc	14,1 cm	4,2 cm	354-5001	RSZ-5001
550 cc	14,4 cm	4,4 cm	354-5501	RSZ-5501
600 cc	14,7 cm	4,5 cm	354-6001	RSZ-6001
700 cc	15,7 cm	4,8 cm	354-7001	RSZ-7001
800 cc	16,6 cm	5,0 cm	354-8001	RSZ-8001

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above

Round Gel Breast Implants, Cohesive I™ CONT'D

SILTEX™ Round Breast Implants, High Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	8,4 cm	3,6 cm	354-4125	RSZ-1254
150 cc	8,9 cm	3,8 cm	354-4150	RSZ-1504
175 cc	9,4 cm	4,0 cm	354-4175	RSZ-1754
200 cc	9,9 cm	4,1 cm	354-4200	RSZ-2004
225 cc	10,2 cm	4,3 cm	354-4225	RSZ-2254
250 cc	10,5 cm	4,5 cm	354-4250	RSZ-2504
275 cc	10,9 cm	4,6 cm	354-4275	RSZ-2754
300 cc	11,1 cm	4,7 cm	354-4300	RSZ-3004
325 cc	11,5 cm	4,8 cm	354-4325	RSZ-3254
350 cc	11,7 cm	4,9 cm	354-4350	RSZ-3504
375 cc	12,0 cm	5,0 cm	354-4375	RSZ-3754
400 cc	12,3 cm	5,1 cm	354-4400	RSZ-4004
425 cc	12,5 cm	5,2 cm	354-4425	RSZ-4254
450 cc	12,7 cm	5,2 cm	354-4450	RSZ-4504
500 cc	13,2 cm	5,4 cm	354-4500	RSZ-5004
550 cc	13,5 cm	5,6 cm	354-4550	RSZ-5504
600 cc	14,0 cm	5,7 cm	354-4600	RSZ-6004
650 cc	14,3 cm	5,8 cm	354-4650	RSZ-6504
700 cc	14,7 cm	6,0 cm	354-4700	RSZ-7004
800 cc	15,4 cm	6,3 cm	354-4800	RSZ-8004

SILTEX™ Round Breast Implants, Ultra High Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
135 cc	8,0 cm	4,3 cm	354-5135	RSZ-5135
160 cc	8,4 cm	4,5 cm	354-5160	RSZ-5160
185 cc	8,6 cm	4,6 cm	354-5185	RSZ-5185
215 cc	8,9 cm	4,7 cm	354-5215	RSZ-5215
240 cc	9,2 cm	4,9 cm	354-5240	RSZ-5240
270 cc	9,5 cm	5,0 cm	354-5270	RSZ-5270
295 cc	9,9 cm	5,2 cm	354-5295	RSZ-5295
320 cc	10,0 cm	5,3 cm	354-5320	RSZ-5320
350 cc	10,3 cm	5,4 cm	354-5350	RSZ-5350
375 cc	10,5 cm	5,5 cm	354-5375	RSZ-5375
400 cc	10,8 cm	5,6 cm	354-5400	RSZ-5400
430 cc	11,1 cm	5,8 cm	354-5430	RSZ-5430
455 cc	11,3 cm	5,9 cm	354-5455	RSZ-5455
480 cc	11,6 cm	6,0 cm	354-5480	RSZ-5480
535 cc	12,1 cm	6,3 cm	354-5535	RSZ-5535
590 cc	12,6 cm	6,5 cm	354-5590	RSZ-5590

*Resterilizable up to 10 times.

MEMORYGEL™ BREAST IMPLANTS

Round Gel Breast Implants, Cohesive II™

SILTEX™ Round Breast Implants, Moderate Profile, Cohesive II™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
100 cc	9,1 cm	2,6 cm	354-1009	RSZ-7100
125 cc	9,3 cm	2,8 cm	354-1259	RSZ-7125
150 cc	10,1 cm	2,4 cm	354-1509	RSZ-7150
175 cc	10,5 cm	2,9 cm	354-1759	RSZ-7175
200 cc	11,2 cm	2,8 cm	354-2009	RSZ-7200
225 cc	11,5 cm	3,0 cm	354-2259	RSZ-7225
250 cc	11,6 cm	3,5 cm	354-2509	RSZ-7250
275 cc	12,5 cm	3,1 cm	354-2759	RSZ-7275
300 cc	12,6 cm	3,6 cm	354-3009	RSZ-7300
325 cc	13,1 cm	3,6 cm	354-3259	RSZ-7325
350 cc	13,4 cm	3,5 cm	354-3509	RSZ-7350
375 cc	13,4 cm	3,9 cm	354-3759	RSZ-7375
400 cc	13,7 cm	3,8 cm	354-4009	RSZ-7400
450 cc	13,7 cm	4,6 cm	354-4509	RSZ-7450
500 cc	14,2 cm	4,8 cm	354-5009	RSZ-7500
550 cc	14,5 cm	4,5 cm	354-5509	RSZ-7550
600 cc	15,3 cm	4,6 cm	354-6009	RSZ-7600
700 cc	16,3 cm	5,0 cm	354-7009	RSZ-7700
800 cc	16,7 cm	5,4 cm	354-8009	RSZ-7800

SILTEX™ Round Breast Implants, Moderate Plus Profile, Cohesive II™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	8,8 cm	2,9 cm	324-5125	RSZ-1251
150 cc	9,4 cm	3,0 cm	324-5150	RSZ-1501
175 cc	10,0 cm	3,2 cm	324-5175	RSZ-1751
200 cc	10,5 cm	3,3 cm	324-5200	RSZ-2001
225 cc	10,9 cm	3,5 cm	324-5225	RSZ-2251
250 cc	11,3 cm	3,6 cm	324-5250	RSZ-2501
275 cc	11,7 cm	3,7 cm	324-5275	RSZ-2751
300 cc	12,0 cm	3,7 cm	324-5300	RSZ-3001
325 cc	12,3 cm	3,8 cm	324-5325	RSZ-3251
350 cc	12,6 cm	3,8 cm	324-5350	RSZ-3501
375 cc	12,9 cm	3,9 cm	324-5375	RSZ-3751
400 cc	13,2 cm	4,0 cm	324-5400	RSZ-4001
450 cc	13,7 cm	4,1 cm	324-5450	RSZ-4501
500 cc	14,1 cm	4,2 cm	324-5500	RSZ-5001
550 cc	14,4 cm	4,4 cm	324-5550	RSZ-5501
600 cc	14,7 cm	4,5 cm	324-5600	RSZ-6001
700 cc	15,7 cm	4,8 cm	324-5700	RSZ-7001
800 cc	16,6 cm	5,0 cm	324-5800	RSZ-8001

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

Round Gel Breast Implants, Cohesive II™ CONT'D

SILTEX™ Round Breast Implants, High Profile, Cohesive II™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	8,4 cm	3,6 cm	324-4125	RSZ-1254
150 cc	8,9 cm	3,8 cm	324-4150	RSZ-1504
175 cc	9,4 cm	4,0 cm	324-4175	RSZ-1754
200 cc	9,9 cm	4,1 cm	324-4200	RSZ-2004
225 cc	10,2 cm	4,3 cm	324-4225	RSZ-2254
250 cc	10,5 cm	4,5 cm	324-4250	RSZ-2504
275 cc	10,9 cm	4,6 cm	324-4275	RSZ-2754
300 cc	11,1 cm	4,7 cm	324-4300	RSZ-3004
325 cc	11,5 cm	4,8 cm	324-4325	RSZ-3254
350 cc	11,7 cm	4,9 cm	324-4350	RSZ-3504
375 cc	12,0 cm	5,0 cm	324-4375	RSZ-3754
400 cc	12,3 cm	5,1 cm	324-4400	RSZ-4004
425 cc	12,5 cm	5,2 cm	324-4425	RSZ-4254
450 cc	12,7 cm	5,2 cm	324-4450	RSZ-4504
500 cc	13,2 cm	5,4 cm	324-4500	RSZ-5004
550 cc	13,5 cm	5,6 cm	324-4550	RSZ-5504
600 cc	14,0 cm	5,7 cm	324-4600	RSZ-6004
650 cc	14,3 cm	5,8 cm	324-4650	RSZ-6504
700 cc	14,7 cm	6,0 cm	324-4700	RSZ-7004
800 cc	15,4 cm	6,3 cm	324-4800	RSZ-8004

*Resterilizable up to 10 times.

CPG™ Gel Breast Implants, Cohesive III™

CPG™ BREAST IMPLANT NAMING LEGEND		
Cohesive Level	Implant Height	Implant Projection
III 3	Tall 3	High 3
	Medium 2	Moderate Plus 2
	Low 1	Moderate 1
KEY: DIGIT 1 = COHESIVE LEVEL DIGIT 2 = HEIGHT OF IMPLANT DIGIT 3 = PROJECTION OF IMPLANT		

UNDERSTANDING THE LOWER POLE ARC

The diagram illustrates the measurement of the Lower Pole Arc (LPA) on a breast implant. It shows a cross-section of the implant with a vertical line on the left representing the height. Brackets indicate the top 50% and bottom 50% of the height, with the midpoint labeled 'MIDWAY POINT'. The bottom of the implant is labeled 'B' and '0,5 cm'. The curved line along the bottom is labeled 'LOWER POLE ARC'. The right side of the implant is labeled 'A'. A note states: 'MEASUREMENT ASSUMES 0,5 CM OF TISSUE THICKNESS'.

LOWER POLE ARC (LPA) =

Length measured along the lower pole of the implant running from the bottom of the implant (Point B) upto the midway point (Point A) including a 0,5 cm of tissue thickness.

Diagrams not to scale, for reference only.

CPG™ Gel Breast Implants, Cohesive III™ CONT'D

CPG™ Implants 331,
Cohesive III™, Tall Height, Moderate Projection

Volume	Width	Height	Projection	Lower Pole Arc A-B	Catalog #	Resterilizable Gel Sizer*
125 cc	9,0 cm	9,2 cm	3,2 cm	7,3 cm	334-0903	RSZ-0903
150 cc	9,5 cm	9,7 cm	3,3 cm	7,7 cm	334-0953	RSZ-0953
175 cc	10,0 cm	10,2 cm	3,4 cm	8,1 cm	334-1003	RSZ-1003
200 cc	10,5 cm	10,7 cm	3,5 cm	8,4 cm	334-1053	RSZ-1053
230 cc	11,0 cm	11,3 cm	3,6 cm	8,8 cm	334-1103	RSZ-1103
265 cc	11,5 cm	11,8 cm	3,7 cm	9,1 cm	334-1153	RSZ-1153
300 cc	12,0 cm	12,3 cm	3,9 cm	9,5 cm	334-1203	RSZ-1203
340 cc	12,5 cm	12,8 cm	4,0 cm	9,9 cm	334-1253	RSZ-1253
380 cc	13,0 cm	13,3 cm	4,1 cm	10,2 cm	334-1303	RSZ-1303
425 cc	13,5 cm	13,8 cm	4,3 cm	10,6 cm	334-1353	RSZ-1353
475 cc	14,0 cm	14,3 cm	4,5 cm	11,0 cm	334-1403	RSZ-1403
530 cc	14,5 cm	14,8 cm	4,6 cm	11,3 cm	334-1453	RSZ-1453
585 cc	15,0 cm	15,3 cm	4,8 cm	11,7 cm	334-1503	RSZ-1503
645 cc	15,5 cm	15,9 cm	5,0 cm	12,1 cm	334-1553	RSZ-1553

CPG™ Implants 321,
Cohesive III™, Medium Height, Moderate Projection

Volume	Width	Height	Projection	Lower Pole Arc A-B	Catalog #	Resterilizable Gel Sizer*
120 cc	9,0 cm	8,5 cm	3,3 cm	7,0 cm	354-0908	RSZ-0908
135 cc	9,5 cm	8,9 cm	3,5 cm	7,3 cm	354-0958	RSZ-0958
155 cc	10,0 cm	9,4 cm	3,7 cm	7,7 cm	354-1008	RSZ-1008
180 cc	10,5 cm	9,9 cm	3,8 cm	8,0 cm	354-1058	RSZ-1058
215 cc	11,0 cm	10,3 cm	3,9 cm	8,2 cm	354-1108	RSZ-1108
245 cc	11,5 cm	10,8 cm	4,0 cm	8,6 cm	354-1158	RSZ-1158
280 cc	12,0 cm	11,3 cm	4,2 cm	8,9 cm	354-1208	RSZ-1208
315 cc	12,5 cm	11,8 cm	4,4 cm	9,2 cm	354-1258	RSZ-1258
355 cc	13,0 cm	12,2 cm	4,6 cm	9,6 cm	354-1308	RSZ-1308
395 cc	13,5 cm	12,7 cm	4,7 cm	9,9 cm	354-1358	RSZ-1358
440 cc	14,0 cm	13,2 cm	4,9 cm	10,3 cm	354-1408	RSZ-1408
480 cc	14,5 cm	13,6 cm	5,0 cm	10,5 cm	354-1458	RSZ-1458
530 cc	15,0 cm	14,1 cm	5,2 cm	10,9 cm	354-1508	RSZ-1508
640 cc	16,0 cm	15,0 cm	5,6 cm	11,5 cm	354-1608	RSZ-1608
775 cc	17,0 cm	16,0 cm	5,9 cm	12,2 cm	354-1708	RSZ-1708

CPG™ Implants 311,
Cohesive III™, Low Height, Moderate Projection

Volume	Width	Height	Projection	Lower Pole Arc A-B	Catalog #	Resterilizable Gel Sizer*
120 cc	9,5 cm	8,5 cm	3,2 cm	6,9 cm	334-0951	RSZ-0095
140 cc	10,0 cm	8,9 cm	3,4 cm	7,2 cm	334-1001	RSZ-0100
160 cc	10,5 cm	9,4 cm	3,6 cm	7,5 cm	334-1051	RSZ-0105
180 cc	11,0 cm	9,8 cm	3,7 cm	7,9 cm	334-1101	RSZ-0110
210 cc	11,5 cm	10,3 cm	3,9 cm	8,2 cm	334-1151	RSZ-0115
235 cc	12,0 cm	10,7 cm	4,1 cm	8,5 cm	334-1201	RSZ-0120
270 cc	12,5 cm	11,2 cm	4,2 cm	8,8 cm	334-1251	RSZ-0125
300 cc	13,0 cm	11,6 cm	4,4 cm	9,2 cm	334-1301	RSZ-0130
335 cc	13,5 cm	12,0 cm	4,6 cm	9,5 cm	334-1351	RSZ-0135
375 cc	14,0 cm	12,5 cm	4,7 cm	9,8 cm	334-1401	RSZ-0140
415 cc	14,5 cm	12,9 cm	4,9 cm	10,1 cm	334-1451	RSZ-0145
460 cc	15,0 cm	13,4 cm	5,1 cm	10,5 cm	334-1501	RSZ-0150
510 cc	15,5 cm	13,8 cm	5,2 cm	10,8 cm	334-1551	RSZ-0155
560 cc	16,0 cm	14,3 cm	5,4 cm	11,1 cm	334-1601	RSZ-0160
615 cc	16,5 cm	14,7 cm	5,6 cm	11,4 cm	334-1651	RSZ-0165

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

MEMORYGEL™ BREAST IMPLANTS

CPG™ Gel Breast Implants, Cohesive III™ CONT'D

CPG™ Implants 332, Cohesive III™, Tall Height, Moderate Plus Projection



Volume	Width	Height	Projection	Lower Pole		Catalog #	Resterilizable Gel Sizer*
				Arc A-B			
145 cc	9,0 cm	9,4 cm	3,8 cm	8,1 cm		334-0909	RSZ-0909
175 cc	9,5 cm	9,9 cm	4,0 cm	8,5 cm		334-0959	RSZ-0959
205 cc	10,0 cm	10,4 cm	4,2 cm	8,9 cm		334-1009	RSZ-1009
235 cc	10,5 cm	10,9 cm	4,4 cm	9,3 cm		334-1059	RSZ-1059
270 cc	11,0 cm	11,5 cm	4,7 cm	9,7 cm		334-1109	RSZ-1109
305 cc	11,5 cm	12,0 cm	4,9 cm	10,1 cm		334-1159	RSZ-1159
350 cc	12,0 cm	12,5 cm	5,1 cm	10,5 cm		334-1209	RSZ-1209
395 cc	12,5 cm	13,0 cm	5,3 cm	10,9 cm		334-1259	RSZ-1259
445 cc	13,0 cm	13,5 cm	5,5 cm	11,3 cm		334-1309	RSZ-1309
495 cc	13,5 cm	14,1 cm	5,7 cm	11,7 cm		334-1359	RSZ-1359
555 cc	14,0 cm	14,6 cm	5,9 cm	12,1 cm		334-1409	RSZ-1409
615 cc	14,5 cm	15,1 cm	6,1 cm	12,6 cm		334-1459	RSZ-1459
680 cc	15,0 cm	15,6 cm	6,3 cm	13,0 cm		334-1509	RSZ-1509

CPG™ Implants 322, Cohesive III™, Medium Height, Moderate Plus Projection



Volume	Width	Height	Projection	Lower Pole		Resterilizable Gel Sizer*
				Arc A-B	Catalog #	
140 cc	9,0 cm	8,5 cm	3,8 cm	7,6 cm	334-0905	RSZ-0905
165 cc	9,5 cm	8,9 cm	4,0 cm	8,0 cm	334-0955	RSZ-0955
195 cc	10,0 cm	9,4 cm	4,2 cm	8,4 cm	334-1005	RSZ-1005
225 cc	10,5 cm	9,9 cm	4,4 cm	8,8 cm	334-1055	RSZ-1055
255 cc	11,0 cm	10,3 cm	4,7 cm	9,2 cm	334-1105	RSZ-1105
295 cc	11,5 cm	10,8 cm	4,9 cm	9,5 cm	334-1155	RSZ-1155
330 cc	12,0 cm	11,3 cm	5,1 cm	9,9 cm	334-1205	RSZ-1205
375 cc	12,5 cm	11,8 cm	5,3 cm	10,3 cm	334-1255	RSZ-1255
420 cc	13,0 cm	12,2 cm	5,5 cm	10,7 cm	334-1305	RSZ-1305
475 cc	13,5 cm	12,7 cm	5,7 cm	11,1 cm	334-1355	RSZ-1355
525 cc	14,0 cm	13,2 cm	5,9 cm	11,4 cm	334-1405	RSZ-1405
585 cc	14,5 cm	13,6 cm	6,1 cm	11,8 cm	334-1455	RSZ-1455
650 cc	15,0 cm	14,1 cm	6,3 cm	12,2 cm	334-1505	RSZ-1505

CPG™ Implants 312, Cohesive III™, Low Height, Moderate Plus Projection



Volume	Width	Height	Projection	Lower Pole		Resterilizable Gel Sizer*
				Arc A-B	Catalog #	
125 cc	9,0 cm	8,0 cm	3,8 cm	7,2 cm	334-0907	RSZ-0907
145 cc	9,5 cm	8,4 cm	4,0 cm	7,5 cm	334-0957	RSZ-0957
170 cc	10,0 cm	8,8 cm	4,2 cm	7,9 cm	334-1007	RSZ-1007
195 cc	10,5 cm	9,3 cm	4,4 cm	8,3 cm	334-1057	RSZ-1057
225 cc	11,0 cm	9,7 cm	4,7 cm	8,6 cm	334-1107	RSZ-1107
255 cc	11,5 cm	10,2 cm	4,9 cm	9,0 cm	334-1157	RSZ-1157
290 cc	12,0 cm	10,6 cm	5,1 cm	9,3 cm	334-1207	RSZ-1207
330 cc	12,5 cm	11,1 cm	5,3 cm	9,7 cm	334-1257	RSZ-1257
370 cc	13,0 cm	11,5 cm	5,5 cm	10,0 cm	334-1307	RSZ-1307
415 cc	13,5 cm	11,9 cm	5,7 cm	10,4 cm	334-1357	RSZ-1357
465 cc	14,0 cm	12,4 cm	5,9 cm	10,8 cm	334-1407	RSZ-1407
515 cc	14,5 cm	12,8 cm	6,1 cm	11,1 cm	334-1457	RSZ-1457
570 cc	15,0 cm	13,3 cm	6,3 cm	11,5 cm	334-1507	RSZ-1507
690 cc	16,0 cm	14,2 cm	6,8 cm	12,2 cm	334-1607	RSZ-1607

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

CPG™ Gel Breast Implants, Cohesive III™ CONT'D

CPG™ Implants 333,
Cohesive III™, Tall Height, High Projection

Volume	Width	Height	Projection	Lower Pole		Resterilizable Gel Sizer*
				Arc A-B	Catalog #	
180 cc	9,0 cm	9,4 cm	4,6 cm	8,4 cm	334-0904	RSZ-4090
215 cc	9,5 cm	9,9 cm	4,8 cm	8,8 cm	334-0954	RSZ-4095
250 cc	10,0 cm	10,4 cm	5,1 cm	9,3 cm	334-1004	RSZ-4100
290 cc	10,5 cm	10,9 cm	5,3 cm	9,7 cm	334-1054	RSZ-4105
330 cc	11,0 cm	11,5 cm	5,6 cm	10,1 cm	334-1104	RSZ-4110
380 cc	11,5 cm	12,0 cm	5,8 cm	10,5 cm	334-1154	RSZ-4115
430 cc	12,0 cm	12,5 cm	6,0 cm	11,0 cm	334-1204	RSZ-4120
485 cc	12,5 cm	13,0 cm	6,2 cm	11,4 cm	334-1254	RSZ-4125
545 cc	13,0 cm	13,5 cm	6,5 cm	11,8 cm	334-1304	RSZ-4130
610 cc	13,5 cm	14,1 cm	6,7 cm	12,2 cm	334-1354	RSZ-4135
680 cc	14,0 cm	14,6 cm	6,9 cm	12,7 cm	334-1404	RSZ-4140
755 cc	14,5 cm	15,1 cm	7,1 cm	13,1 cm	334-1454	RSZ-4145

CPG™ Implants 323,
Cohesive III™, Medium Height, High Projection

Volume	Width	Height	Projection	Lower Pole		Resterilizable Gel Sizer*
				Arc A-B	Catalog #	
165 cc	9,0 cm	8,5 cm	4,6 cm	8,1 cm	334-0902	RSZ-0902
195 cc	9,5 cm	8,9 cm	4,8 cm	8,5 cm	334-0952	RSZ-0952
225 cc	10,0 cm	9,4 cm	5,1 cm	8,9 cm	334-1002	RSZ-1002
260 cc	10,5 cm	9,9 cm	5,3 cm	9,3 cm	334-1052	RSZ-1052
300 cc	11,0 cm	10,3 cm	5,6 cm	9,7 cm	334-1102	RSZ-1102
345 cc	11,5 cm	10,8 cm	5,8 cm	10,1 cm	334-1152	RSZ-1152
390 cc	12,0 cm	11,3 cm	6,0 cm	10,5 cm	334-1202	RSZ-1202
440 cc	12,5 cm	11,8 cm	6,2 cm	10,9 cm	334-1252	RSZ-1252
495 cc	13,0 cm	12,2 cm	6,5 cm	11,3 cm	334-1302	RSZ-1302
555 cc	13,5 cm	12,7 cm	6,7 cm	11,8 cm	334-1352	RSZ-1352
620 cc	14,0 cm	13,2 cm	6,9 cm	12,2 cm	334-1402	RSZ-1402
685 cc	14,5 cm	13,6 cm	7,1 cm	12,6 cm	334-1452	RSZ-1452

CPG™ Implants 313,
Cohesive III™, Low Height, High Projection

Volume	Width	Height	Projection	Lower Pole		Resterilizable Gel Sizer*
				Arc A-B	Catalog #	
130 cc	9,0 cm	8,1 cm	4,4 cm	7,6 cm	334-0906	RSZ-0906
155 cc	9,5 cm	8,6 cm	4,5 cm	7,8 cm	334-0956	RSZ-0956
180 cc	10,0 cm	9,0 cm	4,6 cm	8,3 cm	334-1006	RSZ-1006
210 cc	10,5 cm	9,5 cm	4,8 cm	8,7 cm	334-1056	RSZ-1056
240 cc	11,0 cm	9,9 cm	4,9 cm	9,1 cm	334-1106	RSZ-1106
270 cc	11,5 cm	10,3 cm	5,0 cm	9,5 cm	334-1156	RSZ-1156
310 cc	12,0 cm	10,8 cm	5,2 cm	9,9 cm	334-1206	RSZ-1206
350 cc	12,5 cm	11,2 cm	5,4 cm	10,2 cm	334-1256	RSZ-1256
395 cc	13,0 cm	11,7 cm	5,6 cm	10,6 cm	334-1306	RSZ-1306
440 cc	13,5 cm	12,2 cm	5,8 cm	11,0 cm	334-1356	RSZ-1356
490 cc	14,0 cm	12,6 cm	6,1 cm	11,4 cm	334-1406	RSZ-1406
545 cc	14,5 cm	13,1 cm	6,3 cm	11,7 cm	334-1456	RSZ-1456
605 cc	15,0 cm	13,5 cm	6,6 cm	12,1 cm	334-1506	RSZ-1506

*Resterilizable up to 10 times.

MEMORYGEL™ BREAST IMPLANTS

MENTOR® CONTOUR PROFILE™ Breast Tissue Expanders / CPG™ Gel Breast Implants Sizing Guide – Two-Stage Breast Reconstruction

STAGE 1

CPX™4 Breast Tissue Expander, and CPX™ 4 Breast Tissue Expander with Suture Tabs

CPX™4 Expander, 8100 Style & CPX™4 Expander, 9100 Style with Suture Tabs – Low Height

Volume	Width	Height	Proj.	Catalog #
250 cc	11,4 cm	8,1 cm	6,1 cm	354-8121 354-9121
350 cc	12,7 cm	9,4 cm	6,5 cm	354-8122 354-9122
450 cc	14,0 cm	10,2 cm	7,1 cm	354-8123 354-9123
550 cc	15,0 cm	10,9 cm	7,4 cm	354-8124 354-9124
650 cc	15,7 cm	11,2 cm	7,9 cm	354-8125 354-9125
750 cc	16,5 cm	11,9 cm	8,1 cm	354-8126 354-9126

CPX™4 Expander, 8200 Style & CPX™4 Expander, 9200 Style with Suture Tabs – Medium Height

Volume	Width	Height	Proj.	Catalog #
275 cc	10,7 cm	9,3 cm	6,2 cm	354-8221 354-9221
350 cc	11,7 cm	10,0 cm	6,6 cm	354-8222 354-9222
450 cc	12,7 cm	10,8 cm	7,0 cm	354-8223 354-9223
550 cc	13,5 cm	11,7 cm	7,4 cm	354-8224 354-9224
650 cc	14,6 cm	12,6 cm	7,6 cm	354-8225 354-9225
800 cc	15,6 cm	13,3 cm	8,0 cm	354-8226 354-9226

CPX™4 Expander, 8300 Style & CPX™4 Expander, 9300 Style with Suture Tabs – Tall Height

Volume	Width	Height	Proj.	Catalog #
250 cc	10,1 cm	10,7 cm	5,6 cm	354-8321 354-9321
350 cc	11,3 cm	11,8 cm	6,0 cm	354-8322 354-9322
450 cc	12,3 cm	12,9 cm	6,5 cm	354-8323 354-9323
550 cc	13,2 cm	13,8 cm	6,9 cm	354-8324 354-9324
650 cc	14,0 cm	14,6 cm	7,3 cm	354-8325 354-9325
750 cc	14,6 cm	15,3 cm	7,6 cm	354-8326 354-9326
850 cc	15,4 cm	15,9 cm	7,9 cm	354-8327 354-9327

Selection of CPG™ Gel Breast

CPG™ Implant 311

Cohesive III™, Low Height, Moderate Projection

Volume	Width	Height	Proj.	Catalog #
180 cc	11,0 cm	9,8 cm	3,7 cm	334-1101
210 cc	11,5 cm	10,3 cm	3,9 cm	334-1151
270 cc	12,5 cm	11,2 cm	4,2 cm	334-1251
300 cc	13,0 cm	11,6 cm	4,4 cm	334-1301
375 cc	14,0 cm	12,5 cm	4,7 cm	334-1401
460 cc	15,0 cm	13,4 cm	5,1 cm	334-1501
510 cc	15,5 cm	13,8 cm	5,2 cm	334-1551
560 cc	16,0 cm	14,3 cm	5,4 cm	334-1601
615 cc	16,5 cm	14,7 cm	5,6 cm	334-1651

CPG™ Implant 311

Cohesive III™, Low Height, Moderate Projection

Volume	Width	Height	Proj.	Catalog #
160 cc	10,5 cm	9,4 cm	3,6 cm	334-1051
180 cc	11,0 cm	9,8 cm	3,7 cm	334-1101
210 cc	11,5 cm	10,3 cm	3,9 cm	334-1151
235 cc	12,0 cm	10,7 cm	4,1 cm	334-1201
270 cc	12,5 cm	11,2 cm	4,2 cm	334-1251
300 cc	13,0 cm	11,6 cm	4,4 cm	334-1301
335 cc	13,5 cm	12,0 cm	4,6 cm	334-1351
415 cc	14,5 cm	12,9 cm	4,9 cm	334-1451
460 cc	15,0 cm	13,4 cm	5,1 cm	334-1501
510 cc	15,5 cm	13,8 cm	5,2 cm	334-1551
560 cc	16,0 cm	14,3 cm	5,4 cm	334-1601

CPG™ Implant 321

Cohesive III™, Medium Height, Moderate Projection

Volume	Width	Height	Proj.	Catalog #
180 cc	10,5 cm	9,9 cm	3,8 cm	354-1058
215 cc	11,0 cm	10,3 cm	3,9 cm	354-1108
245 cc	11,5 cm	10,8 cm	4,0 cm	354-1158
280 cc	12,0 cm	11,3 cm	4,2 cm	354-1208
315 cc	12,5 cm	11,8 cm	4,4 cm	354-1258
355 cc	13,0 cm	12,2 cm	4,6 cm	354-1308
395 cc	13,5 cm	12,7 cm	4,7 cm	354-1358
480 cc	14,5 cm	13,6 cm	5,0 cm	354-1458
530 cc	15,0 cm	14,1 cm	5,2 cm	354-1508
640 cc	16,0 cm	15,0 cm	5,6 cm	354-1608

CPG™ Implant 331

Cohesive III™, Tall Height, Moderate Projection

Volume	Width	Height	Proj.	Catalog #
175 cc	10,0 cm	10,2 cm	3,4 cm	334-1003
200 cc	10,5 cm	10,7 cm	3,5 cm	334-1053
230 cc	11,0 cm	11,3 cm	3,6 cm	334-1103
265 cc	11,5 cm	11,8 cm	3,7 cm	334-1153
300 cc	12,0 cm	12,3 cm	3,9 cm	334-1203
340 cc	12,5 cm	12,8 cm	4,0 cm	334-1253
380 cc	13,0 cm	13,3 cm	4,1 cm	334-1303
425 cc	13,5 cm	13,8 cm	4,3 cm	334-1353
475 cc	14,0 cm	14,3 cm	4,5 cm	334-1403
530 cc	14,5 cm	14,8 cm	4,6 cm	334-1453
585 cc	15,0 cm	15,3 cm	4,8 cm	334-1503
645 cc	15,5 cm	15,9 cm	5,0 cm	334-1553

MEMORYGEL™ BREAST IMPLANTS

*PLEASE NOTE: The CONTOUR PROFILE™ Tissue Expander Injection Dome is not MRI compatible. For detailed indications, contraindications, warning and precautions associated with the use of this device, please refer to the Product Insert Data Sheet. Individual implant dimensions may vary slightly in products of this type. Not all units conform exactly to the dimensions noted.

STAGE 2

Implant to use as a long-term implant after expansion is complete.

CPG™ Implant 312

Cohesive III™, Low Height, Moderate Plus Projection

Volume	Width	Height	Proj.	Catalog #
225 cc	11,0 cm	9,7 cm	4,7 cm	334-1107
255 cc	11,5 cm	10,2 cm	4,9 cm	334-1157
330 cc	12,5 cm	11,1 cm	5,3 cm	334-1257
370 cc	13,0 cm	11,5 cm	5,5 cm	334-1307
465 cc	14,0 cm	12,4 cm	5,9 cm	334-1407
570 cc	15,0 cm	13,3 cm	6,3 cm	334-1507
690 cc	16,0 cm	14,2 cm	6,8 cm	334-1607

CPG™ Implant 312

Cohesive III™, Low Height, Moderate Plus Projection

Volume	Width	Height	Proj.	Catalog #
195 cc	10,5 cm	9,3 cm	4,4 cm	334-1057
225 cc	11,0 cm	9,7 cm	4,7 cm	334-1107
255 cc	11,5 cm	10,2 cm	4,9 cm	334-1157
290 cc	12,0 cm	10,6 cm	5,1 cm	334-1207
330 cc	12,5 cm	11,1 cm	5,3 cm	334-1257
370 cc	13,0 cm	11,5 cm	5,5 cm	334-1307
415 cc	13,5 cm	11,9 cm	5,7 cm	334-1357
515 cc	14,5 cm	12,8 cm	6,1 cm	334-1457
570 cc	15,0 cm	13,3 cm	6,3 cm	334-1507
690 cc	16,0 cm	14,2 cm	6,8 cm	334-1607

CPG™ Implant 322

Cohesive III™, Medium Height, Moderate Plus Projection

Volume	Width	Height	Proj.	Catalog #
225 cc	10,5 cm	9,9 cm	4,4 cm	334-1055
255 cc	11,0 cm	10,3 cm	4,7 cm	334-1105
295 cc	11,5 cm	10,8 cm	4,9 cm	334-1155
330 cc	12,0 cm	11,3 cm	5,1 cm	334-1205
375 cc	12,5 cm	11,8 cm	5,3 cm	334-1255
420 cc	13,0 cm	12,2 cm	5,5 cm	334-1305
475 cc	13,5 cm	12,7 cm	5,7 cm	334-1355
585 cc	14,5 cm	13,6 cm	6,1 cm	334-1455
650 cc	15,0 cm	14,1 cm	6,3 cm	334-1505

CPG™ Implant 332

Cohesive III™, Tall Height, Moderate Plus Projection

Volume	Width	Height	Proj.	Catalog #
205 cc	10,0 cm	10,4 cm	4,2 cm	334-1009
270 cc	11,0 cm	11,5 cm	4,7 cm	334-1109
305 cc	11,5 cm	12,0 cm	4,9 cm	334-1159
350 cc	12,0 cm	12,5 cm	5,1 cm	334-1209
395 cc	12,5 cm	13,0 cm	5,3 cm	334-1259
445 cc	13,0 cm	13,5 cm	5,5 cm	334-1309
495 cc	13,5 cm	14,1 cm	5,7 cm	334-1359
555 cc	14,0 cm	14,6 cm	5,9 cm	334-1409
615 cc	14,5 cm	15,1 cm	6,1 cm	334-1459
680 cc	15,0 cm	15,6 cm	6,3 cm	334-1509

CPG™ Implant 313

Cohesive III™, Low Height, High Projection

Volume	Width	Height	Proj.	Catalog #
240 cc	11,0 cm	9,9 cm	4,9 cm	334-1106
270 cc	11,5 cm	10,3 cm	5,0 cm	334-1156
350 cc	12,5 cm	11,2 cm	5,4 cm	334-1256
395 cc	13,0 cm	11,7 cm	5,6 cm	334-1306
490 cc	14,0 cm	12,6 cm	6,1 cm	334-1406
605 cc	15,0 cm	13,5 cm	6,6 cm	334-1506

CPG™ Implant 313

Cohesive III™, Low Height, High Projection

Volume	Width	Height	Proj.	Catalog #
210 cc	10,5 cm	9,5 cm	4,8 cm	334-1056
240 cc	11,0 cm	9,9 cm	4,9 cm	334-1106
270 cc	11,5 cm	10,3 cm	5,0 cm	334-1156
310 cc	12,0 cm	10,8 cm	5,2 cm	334-1206
350 cc	12,5 cm	11,2 cm	5,4 cm	334-1256
395 cc	13,0 cm	11,7 cm	5,6 cm	334-1306
440 cc	13,5 cm	12,2 cm	5,8 cm	334-1356
545 cc	14,5 cm	13,1 cm	6,3 cm	334-1456
606 cc	15,0 cm	13,5 cm	6,6 cm	334-1506

CPG™ Implant 323

Cohesive III™, Medium Height, High Projection

Volume	Width	Height	Proj.	Catalog #
260 cc	10,5 cm	9,9 cm	5,3 cm	334-1052
300 cc	11,0 cm	10,3 cm	5,6 cm	334-1102
345 cc	11,5 cm	10,8 cm	5,8 cm	334-1152
390 cc	12,0 cm	11,3 cm	6,0 cm	334-1202
440 cc	12,5 cm	11,8 cm	6,2 cm	334-1252
495 cc	13,0 cm	12,2 cm	6,5 cm	334-1302
555 cc	13,5 cm	12,7 cm	6,7 cm	334-1352
685 cc	14,5 cm	13,6 cm	7,1 cm	334-1452

CPG™ Implant 333

Cohesive III™, Tall Height, High Projection


Volume	Width	Height	Proj.	Catalog #
250 cc	10,0 cm	10,4 cm	5,1 cm	334-1004
330 cc	11,0 cm	11,5 cm	5,6 cm	334-1104
380 cc	11,5 cm	12,0 cm	5,8 cm	334-1154
430 cc	12,0 cm	12,5 cm	6,0 cm	334-1204
485 cc	12,5 cm	13,0 cm	6,2 cm	334-1254
545 cc	13,0 cm	13,5 cm	6,5 cm	334-1304
610 cc	13,5 cm	14,1 cm	6,7 cm	334-1354
680 cc	14,0 cm	14,6 cm	6,9 cm	334-1404
755 cc	14,5 cm	15,1 cm	7,1 cm	334-1454

MEMORYGEL™ BREAST IMPLANT SIZERS

BECKER™ Expander/Breast Implants

SILTEX™ Round BECKER™ 25
Expander/Implants, Cohesive I™


(25% Silicone Gel, Cohesive I™, in outer lumen, 75% Saline in inner lumen)



Volume	Gel Volume	TEMPORARY OVEREXPANSION VOLUMES		Total Saline	Total Gel-Saline	Diameter	Proj.	Cat.#
		Max Saline	Total Gel-Saline					
150 cc	40 cc	185 cc	225 cc	85-150 cc	125-190 cc	10,7 cm	2,4 cm	354-1500
200 cc	50 cc	250 cc	300 cc	125-200 cc	175-250 cc	11,6 cm	2,7 cm	354-2000
250 cc	60 cc	315 cc	375 cc	165-255 cc	225-315 cc	12,1 cm	3,0 cm	354-2500
300 cc	75 cc	375 cc	450 cc	200-300 cc	275-375 cc	12,9 cm	3,0 cm	354-3000
350 cc	90 cc	435 cc	525 cc	235-350 cc	325-440 cc	14,1 cm	3,0 cm	354-3500
400 cc	100 cc	500 cc	600 cc	275-400 cc	375-500 cc	13,9 cm	3,4 cm	354-4000
500 cc	125 cc	625 cc	750 cc	350-500 cc	475-625 cc	14,7 cm	3,7 cm	354-5000
600 cc	150 cc	750 cc	900 cc	425-600 cc	575-750 cc	15,6 cm	4,2 cm	354-6000
700 cc	175 cc	875 cc	1050 cc	500-700 cc	675-875 cc	16,8 cm	4,3 cm	354-7000
800 cc	200 cc	1000 cc	1200 cc	575-800 cc	775-1000 cc	17,2 cm	4,6 cm	354-8000

SILTEX™ Round BECKER™ 50
Expander/Implants, Cohesive I™

(50% Silicone Gel, Cohesive I™, in outer lumen, 50% Saline in inner lumen)



Volume	Gel Volume	Total Saline	Total Gel- Saline	Diameter	Proj.	Catalog #
300 cc	150 cc	150-200 cc	300-350 cc	12,9 cm	3,0 cm	354-1515
400 cc	200 cc	200-300 cc	400-500 cc	13,9 cm	3,4 cm	354-2020
500 cc	250 cc	250-350 cc	500-600 cc	14,7 cm	3,7 cm	354-2525
600 cc	300 cc	300-425 cc	600-725 cc	15,6 cm	4,2 cm	354-3030
700 cc	350 cc	350-500 cc	700-850 cc	16,0 cm	5,0 cm	354-3535

SILTEX™ CONTOUR PROFILE™
BECKER™ 35 Expander/Implants, Cohesive II™

(35% Silicone Gel, Cohesive II™, in outer lumen, 65% Saline in inner lumen)



Volume	Gel Volume	TEMPORARY OVEREXPANSION VOLUMES		Total Saline	Total Gel- Saline	Width	Height	Proj.	Cat.#
		Max Saline	Total Gel-Saline						
145 cc	50 cc	120 cc	170 cc	85-95 cc	135-145 cc	9,5 cm	8,3 cm	4,2 cm	324-0955
195 cc	70 cc	155 cc	225 cc	110-125 cc	180-195 cc	10,5 cm	9,1 cm	4,6 cm	324-1055
255 cc	90 cc	205 cc	295 cc	145-165 cc	235-255 cc	11,5 cm	10,0 cm	5,1 cm	324-1155
290 cc	100 cc	235 cc	335 cc	170-190 cc	270-290 cc	12,0 cm	10,5 cm	5,3 cm	324-1205
325 cc	110 cc	270 cc	380 cc	190-215 cc	300-325 cc	12,5 cm	10,9 cm	5,5 cm	324-1255
365 cc	125 cc	300 cc	425 cc	215-240 cc	340-365 cc	13,0 cm	11,3 cm	5,7 cm	324-1305
400 cc	130 cc	335 cc	465 cc	240-270 cc	370-400 cc	13,5 cm	11,8 cm	5,9 cm	324-1355
460 cc	160 cc	375 cc	535 cc	270-300 cc	430-460 cc	14,0 cm	12,2 cm	6,2 cm	324-1405
565 cc	195 cc	460 cc	655 cc	330-370 cc	525-565 cc	15,0 cm	13,1 cm	6,6 cm	324-1505
685 cc	240 cc	555 cc	795 cc	400-445 cc	640-685 cc	16,0 cm	13,9 cm	7,0 cm	324-1605

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

Sterile Resterilizable Gel Breast Implant Sizers*

Round Gel Breast Implant Sizers, Moderate Profile



Volume	Diameter	Projection	Catalog #
125 cc	10,0 cm	2,2 cm	RSZ-7125
150 cc	10,3 cm	2,3 cm	RSZ-7150
175 cc	11,2 cm	2,4 cm	RSZ-7175
200 cc	11,7 cm	2,5 cm	RSZ-7200
225 cc	12,2 cm	2,6 cm	RSZ-7225
250 cc	12,3 cm	2,8 cm	RSZ-7250
275 cc	13,2 cm	2,9 cm	RSZ-7275
300 cc	13,5 cm	3,0 cm	RSZ-7300
325 cc	13,9 cm	3,0 cm	RSZ-7325
350 cc	14,2 cm	3,1 cm	RSZ-7350
375 cc	14,4 cm	3,2 cm	RSZ-7375
400 cc	14,5 cm	3,2 cm	RSZ-7400
450 cc	14,9 cm	3,2 cm	RSZ-7450
500 cc	15,2 cm	3,4 cm	RSZ-7500
550 cc	15,9 cm	3,6 cm	RSZ-7550
600 cc	16,5 cm	3,7 cm	RSZ-7600
700 cc	17,4 cm	3,9 cm	RSZ-7700
800 cc	18,2 cm	4,1 cm	RSZ-7800

Round Gel Breast Implant Sizers, Moderate Plus Profile



Volume	Diameter	Projection	Catalog #
125 cc	8,9 cm	2,8 cm	RSZ-1251
150 cc	9,5 cm	2,9 cm	RSZ-1501
175 cc	10,0 cm	3,1 cm	RSZ-1751
200 cc	10,5 cm	3,2 cm	RSZ-2001
225 cc	10,9 cm	3,3 cm	RSZ-2251
250 cc	11,3 cm	3,4 cm	RSZ-2501
275 cc	11,7 cm	3,5 cm	RSZ-2751
300 cc	12,0 cm	3,6 cm	RSZ-3001
325 cc	12,3 cm	3,8 cm	RSZ-3251
350 cc	12,5 cm	3,9 cm	RSZ-3501
375 cc	12,8 cm	4,0 cm	RSZ-3751
400 cc	13,1 cm	4,0 cm	RSZ-4001
450 cc	13,6 cm	4,2 cm	RSZ-4501
500 cc	14,1 cm	4,3 cm	RSZ-5001
550 cc	14,6 cm	4,5 cm	RSZ-5501
600 cc	15,0 cm	4,6 cm	RSZ-6001
700 cc	15,8 cm	4,9 cm	RSZ-7001
800 cc	16,5 cm	5,1 cm	RSZ-8001

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

MEMORYGEL™ BREAST IMPLANT SIZERS

Sterile Resterilizable Gel Breast Implant Sizers* CONT'D

Round Gel Breast Implant Sizers, High Profile



Volume	Diameter	Projection	Catalog #
125 cc	8,3 cm	3,5 cm	RSZ-1254
150 cc	8,8 cm	3,7 cm	RSZ-1504
175 cc	9,3 cm	3,9 cm	RSZ-1754
200 cc	9,7 cm	4,0 cm	RSZ-2004
225 cc	10,1 cm	4,2 cm	RSZ-2254
250 cc	10,5 cm	4,3 cm	RSZ-2504
275 cc	10,8 cm	4,4 cm	RSZ-2754
300 cc	11,1 cm	4,5 cm	RSZ-3004
325 cc	11,4 cm	4,6 cm	RSZ-3254
350 cc	11,7 cm	4,8 cm	RSZ-3504
375 cc	12,0 cm	4,8 cm	RSZ-3754
400 cc	12,2 cm	5,0 cm	RSZ-4004
425 cc	12,5 cm	5,0 cm	RSZ-4254
450 cc	12,8 cm	5,1 cm	RSZ-4504
500 cc	13,2 cm	5,3 cm	RSZ-5004
550 cc	13,6 cm	5,5 cm	RSZ-5504
600 cc	14,0 cm	5,6 cm	RSZ-6004
650 cc	14,4 cm	5,7 cm	RSZ-6504
700 cc	14,8 cm	5,8 cm	RSZ-7004
800 cc	15,5 cm	6,0 cm	RSZ-8004

Round Gel Breast Implant Sizers, Ultra High Profile



Volume	Diameter	Projection	Catalog #
135 cc	7,8 cm	4,1 cm	RSZ-5135
160 cc	8,2 cm	4,3 cm	RSZ-5160
185 cc	8,4 cm	4,4 cm	RSZ-5185
215 cc	8,7 cm	4,5 cm	RSZ-5215
240 cc	9,0 cm	4,7 cm	RSZ-5240
270 cc	9,3 cm	4,8 cm	RSZ-5270
295 cc	9,7 cm	5,0 cm	RSZ-5295
320 cc	9,8 cm	5,1 cm	RSZ-5320
350 cc	10,1 cm	5,2 cm	RSZ-5350
375 cc	10,4 cm	5,3 cm	RSZ-5375
400 cc	10,6 cm	5,4 cm	RSZ-5400
430 cc	10,9 cm	5,6 cm	RSZ-5430
455 cc	11,2 cm	5,7 cm	RSZ-5455
480 cc	11,4 cm	5,8 cm	RSZ-5480
535 cc	12,0 cm	6,1 cm	RSZ-5535
590 cc	12,5 cm	6,3 cm	RSZ-5590

*Resterilizable up to 10 times.

Sterile Resterilizable Gel Breast Implant Sizers* CONT'D

CPG™ GEL BREAST IMPLANT SIZERS

CPG™ Implant Sizers 331 Cohesive III™, Tall Height, Moderate Projection



Volume	Width	Height	Projection	Lower Pole		Catalog #
				Arc A-B		
125 cc	9,0 cm	9,2 cm	3,2 cm	7,3 cm		RSZ-0903
150 cc	9,5 cm	9,7 cm	3,3 cm	7,7 cm		RSZ-0953
175 cc	10,0 cm	10,2 cm	3,4 cm	8,1 cm		RSZ-1003
200 cc	10,5 cm	10,7 cm	3,5 cm	8,4 cm		RSZ-1053
230 cc	11,0 cm	11,3 cm	3,6 cm	8,8 cm		RSZ-1103
265 cc	11,5 cm	11,8 cm	3,7 cm	9,1 cm		RSZ-1153
300 cc	12,0 cm	12,3 cm	3,9 cm	9,5 cm		RSZ-1203
340 cc	12,5 cm	12,8 cm	4,0 cm	9,9 cm		RSZ-1253
380 cc	13,0 cm	13,3 cm	4,1 cm	10,2 cm		RSZ-1303
425 cc	13,5 cm	13,8 cm	4,3 cm	10,6 cm		RSZ-1353
475 cc	14,0 cm	14,3 cm	4,5 cm	11,0 cm		RSZ-1403
530 cc	14,5 cm	14,8 cm	4,6 cm	11,3 cm		RSZ-1453
585 cc	15,0 cm	15,3 cm	4,0 cm	11,7 cm		RSZ-1503
645 cc	15,5 cm	15,9 cm	5,0 cm	12,1 cm		RSZ-1553

CPG™ Implant Sizers 321 Cohesive III™, Medium Height, Moderate Projection



Volume	Width	Height	Projection	Lower Pole		Catalog #
				Arc A-B		
120 cc	9,0 cm	8,5 cm	3,3 cm	7,0 cm		RSZ-0908
135 cc	9,5 cm	8,9 cm	3,5 cm	7,3 cm		RSZ-0958
155 cc	10,0 cm	9,4 cm	3,7 cm	7,7 cm		RSZ-1008
180 cc	10,5 cm	9,9 cm	3,8 cm	8,0 cm		RSZ-1058
215 cc	11,0 cm	10,3 cm	3,9 cm	8,2 cm		RSZ-1108
245 cc	11,5 cm	10,8 cm	4,0 cm	8,6 cm		RSZ-1158
280 cc	12,0 cm	11,3 cm	4,2 cm	8,9 cm		RSZ-1208
315 cc	12,5 cm	11,8 cm	4,4 cm	9,2 cm		RSZ-1258
355 cc	13,0 cm	12,2 cm	4,6 cm	9,6 cm		RSZ-1308
395 cc	13,5 cm	12,7 cm	4,7 cm	9,9 cm		RSZ-1358
440 cc	14,0 cm	13,2 cm	4,9 cm	10,3 cm		RSZ-1408
480 cc	14,5 cm	13,6 cm	5,0 cm	10,5 cm		RSZ-1458
530 cc	15,0 cm	14,1 cm	5,2 cm	10,9 cm		RSZ-1508
640 cc	16,0 cm	15,0 cm	5,6 cm	11,5 cm		RSZ-1608
775 cc	17,0 cm	16,0 cm	5,9 cm	12,2 cm		RSZ-1708

CPG™ Implant Sizers 311 Cohesive III™, Low Height, Moderate Projection



Volume	Width	Height	Projection	Lower Pole		Catalog #
				Arc A-B		
120 cc	9,5 cm	8,5 cm	3,2 cm	6,9 cm		RSZ-0095
140 cc	10,0 cm	8,9 cm	3,4 cm	7,2 cm		RSZ-0100
160 cc	10,5 cm	9,4 cm	3,6 cm	7,5 cm		RSZ-0105
180 cc	11,0 cm	9,8 cm	3,7 cm	7,9 cm		RSZ-0110
210 cc	11,5 cm	10,3 cm	3,9 cm	8,2 cm		RSZ-0115
235 cc	12,0 cm	10,7 cm	4,1 cm	8,5 cm		RSZ-0120
270 cc	12,5 cm	11,2 cm	4,2 cm	8,8 cm		RSZ-0125
300 cc	13,0 cm	11,6 cm	4,4 cm	9,2 cm		RSZ-0130
335 cc	13,5 cm	12,0 cm	4,6 cm	9,5 cm		RSZ-0135
375 cc	14,0 cm	12,5 cm	4,7 cm	9,8 cm		RSZ-0140
415 cc	14,5 cm	12,9 cm	4,9 cm	10,1 cm		RSZ-0145
460 cc	15,0 cm	13,4 cm	5,1 cm	10,5 cm		RSZ-0150
510 cc	15,5 cm	13,8 cm	5,2 cm	10,8 cm		RSZ-0155
560 cc	16,0 cm	14,3 cm	5,4 cm	11,1 cm		RSZ-0160
615 cc	16,5 cm	14,7 cm	5,6 cm	11,4 cm		RSZ-0165

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above

MEMORYGEL™ BREAST IMPLANT SIZERS

Sterile Resterilizable Gel Breast Implant Sizers* CONT'D

CPG™ GEL BREAST IMPLANT SIZERS

CPG™ Implant Sizers 332

Cohesive III™, Tall Height, Moderate Plus Projection



Volume	Width	Height	Projection	Lower Pole		Catalog #
				Arc A-B		
145 cc	9,0 cm	9,4 cm	3,8 cm	8,1 cm		RSZ-0909
175 cc	9,5 cm	9,9 cm	4,0 cm	8,5 cm		RSZ-0959
205 cc	10,0 cm	10,4 cm	4,2 cm	8,9 cm		RSZ-1009
235 cc	10,5 cm	10,9 cm	4,4 cm	9,3 cm		RSZ-1059
270 cc	11,0 cm	11,5 cm	4,7 cm	9,7 cm		RSZ-1109
305 cc	11,5 cm	12,0 cm	4,9 cm	10,1 cm		RSZ-1159
350 cc	12,0 cm	12,5 cm	5,1 cm	10,5 cm		RSZ-1209
395 cc	12,5 cm	13,0 cm	5,3 cm	10,9 cm		RSZ-1259
445 cc	13,0 cm	13,5 cm	5,5 cm	11,3 cm		RSZ-1309
495 cc	13,5 cm	14,1 cm	5,7 cm	11,7 cm		RSZ-1359
555 cc	14,0 cm	14,6 cm	5,9 cm	12,1 cm		RSZ-1409
615 cc	14,5 cm	15,1 cm	6,1 cm	12,6 cm		RSZ-1459
680 cc	15,0 cm	15,6 cm	6,3 cm	13,0 cm		RSZ-1509

CPG™ Implant Sizers 322

Cohesive III™, Medium Height, Moderate Plus Projection



Volume	Width	Height	Projection	Lower Pole		Catalog #
				Arc A-B		
140 cc	9,0 cm	8,5 cm	3,8 cm	7,6 cm		RSZ-0905
165 cc	9,5 cm	8,9 cm	4,0 cm	8,0 cm		RSZ-0955
195 cc	10,0 cm	9,4 cm	4,2 cm	8,4 cm		RSZ-1005
225 cc	10,5 cm	9,9 cm	4,4 cm	8,8 cm		RSZ-1055
255 cc	11,0 cm	10,3 cm	4,7 cm	9,2 cm		RSZ-1105
295 cc	11,5 cm	10,8 cm	4,9 cm	9,5 cm		RSZ-1155
330 cc	12,0 cm	11,3 cm	5,1 cm	9,9 cm		RSZ-1205
375 cc	12,5 cm	11,8 cm	5,3 cm	10,3 cm		RSZ-1255
420 cc	13,0 cm	12,2 cm	5,5 cm	10,7 cm		RSZ-1305
475 cc	13,5 cm	12,7 cm	5,7 cm	11,1 cm		RSZ-1355
525 cc	14,0 cm	13,2 cm	5,9 cm	11,4 cm		RSZ-1405
585 cc	14,5 cm	13,6 cm	6,1 cm	11,8 cm		RSZ-1455
650 cc	15,0 cm	14,1 cm	6,3 cm	12,2 cm		RSZ-1505

CPG™ Implant Sizers 312

Cohesive III™, Low Height, Moderate Plus Projection



Volume	Width	Height	Projection.	Lower Pole		Catalog #
				Arc A-B		
125 cc	9,0 cm	8,0 cm	3,8 cm	7,2 cm		RSZ-0907
145 cc	9,5 cm	8,4 cm	4,0 cm	7,5 cm		RSZ-0957
170 cc	10,0 cm	8,8 cm	4,2 cm	7,9 cm		RSZ-1007
195 cc	10,5 cm	9,3 cm	4,4 cm	8,3 cm		RSZ-1057
225 cc	11,0 cm	9,7 cm	4,7 cm	8,6 cm		RSZ-1107
255 cc	11,5 cm	10,2 cm	4,9 cm	9,0 cm		RSZ-1157
290 cc	12,0 cm	10,6 cm	5,1 cm	9,3 cm		RSZ-1207
330 cc	12,5 cm	11,1 cm	5,3 cm	9,7 cm		RSZ-1257
370 cc	13,0 cm	11,5 cm	5,5 cm	10,0 cm		RSZ-1307
415 cc	13,5 cm	11,9 cm	5,7 cm	10,4 cm		RSZ-1357
465 cc	14,0 cm	12,4 cm	5,9 cm	10,8 cm		RSZ-1407
515 cc	14,5 cm	12,8 cm	6,1 cm	11,1 cm		RSZ-1457
570 cc	15,0 cm	13,3 cm	6,3 cm	11,5 cm		RSZ-1507
690 cc	16,0 cm	14,2 cm	6,8 cm	12,2 cm		RSZ-1607

*Resterilizable up to 10 times.

Sterile Resterilizable Gel Breast Implant Sizers* CONT'D

CPG™ GEL BREAST IMPLANT SIZERS

CPG™ Implant Sizers 333 Cohesive III™, Tall Height, High Projection



Volume	Width	Height	Projection	Lower Pole		Catalog #
				Arc A-B		
180 cc	9,0 cm	9,4 cm	4,6 cm	8,4 cm		RSZ-4090
215 cc	9,5 cm	9,9 cm	4,8 cm	8,8 cm		RSZ-4095
250 cc	10,0 cm	10,4 cm	5,1 cm	9,3 cm		RSZ-4100
290 cc	10,5 cm	10,9 cm	5,3 cm	9,7 cm		RSZ-4105
330 cc	11,0 cm	11,5 cm	5,6 cm	10,1 cm		RSZ-4110
380 cc	11,5 cm	12,0 cm	5,8 cm	10,5 cm		RSZ-4115
430 cc	12,0 cm	12,5 cm	6,0 cm	11,0 cm		RSZ-4120
485 cc	12,5 cm	13,0 cm	6,2 cm	11,4 cm		RSZ-4125
545 cc	13,0 cm	13,5 cm	6,5 cm	11,8 cm		RSZ-4130
610 cc	13,5 cm	14,1 cm	6,7 cm	12,2 cm		RSZ-4135
680 cc	14,0 cm	14,6 cm	6,9 cm	12,7 cm		RSZ-4140
755 cc	14,5 cm	15,1 cm	7,1 cm	13,1 cm		RSZ-4145

CPG™ Implant Sizers 323 Cohesive III™, Medium Height, High Projection



Volume	Width	Height	Projection	Lower Pole		Catalog #
				Arc A-B		
165 cc	9,0 cm	8,5 cm	4,6 cm	8,1 cm		RSZ-0902
195 cc	9,5 cm	8,9 cm	4,8 cm	8,5 cm		RSZ-0952
225 cc	10,0 cm	9,4 cm	5,1 cm	8,9 cm		RSZ-1002
260 cc	10,5 cm	9,9 cm	5,3 cm	9,3 cm		RSZ-1052
300 cc	11,0 cm	10,3 cm	5,6 cm	9,7 cm		RSZ-1102
345 cc	11,5 cm	10,8 cm	5,8 cm	10,1 cm		RSZ-1152
390 cc	12,0 cm	11,3 cm	6,0 cm	10,5 cm		RSZ-1202
440 cc	12,5 cm	11,8 cm	6,2 cm	10,9 cm		RSZ-1252
495 cc	13,0 cm	12,2 cm	6,5 cm	11,3 cm		RSZ-1302
555 cc	13,5 cm	12,7 cm	6,7 cm	11,8 cm		RSZ-1352
620 cc	14,0 cm	13,2 cm	6,9 cm	12,2 cm		RSZ-1402
685 cc	14,5 cm	13,6 cm	7,1 cm	12,6 cm		RSZ-1452

CPG™ Implant Sizers 313 Cohesive III™, Low Height, High Projection



Volume	Width	Height	Projection	Lower Pole		Catalog #
				Arc A-B		
130 cc	9,0 cm	8,1 cm	4,4 cm	7,6 cm		RSZ-0906
155 cc	9,5 cm	8,6 cm	4,5 cm	7,8 cm		RSZ-0956
180 cc	10,0 cm	9,0 cm	4,6 cm	8,3 cm		RSZ-1006
210 cc	10,5 cm	9,5 cm	4,8 cm	8,7 cm		RSZ-1056
240 cc	11,0 cm	9,9 cm	4,9 cm	9,1 cm		RSZ-1106
270 cc	11,5 cm	10,3 cm	5,0 cm	9,5 cm		RSZ-1156
310 cc	12,0 cm	10,8 cm	5,2 cm	9,9 cm		RSZ-1206
350 cc	12,5 cm	11,2 cm	5,4 cm	10,2 cm		RSZ-1256
395 cc	13,0 cm	11,7 cm	5,6 cm	10,6 cm		RSZ-1306
440 cc	13,5 cm	12,2 cm	5,8 cm	11,0 cm		RSZ-1356
490 cc	14,0 cm	12,6 cm	6,1 cm	11,4 cm		RSZ-1406
545 cc	14,5 cm	13,1 cm	6,3 cm	11,7 cm		RSZ-1456
605 cc	15,0 cm	13,5 cm	6,6 cm	12,1 cm		RSZ-1506

*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

TISSUE EXPANDERS

CONTOUR PROFILE™ Tissue Expanders

CPX™ 4 Breast Tissue Expander, Low Height, Style 8100*



Volume	Width	Height	Projection	Catalog #
250 cc	11,4 cm	8,1 cm	6,1 cm	354-8121
350 cc	12,7 cm	9,4 cm	6,5 cm	354-8122
450 cc	14,0 cm	10,2 cm	7,1 cm	354-8123
550 cc	15,0 cm	10,9 cm	7,4 cm	354-8124
650 cc	15,7 cm	11,2 cm	7,9 cm	354-8125
750 cc	16,5 cm	11,9 cm	8,1 cm	354-8126

CPX™ 4 Breast Tissue Expander, Medium Height, Style 8200*



Volume	Width	Height	Projection	Catalog #
275 cc	10,7 cm	9,3 cm	6,2 cm	354-8221
350 cc	11,7 cm	10,0 cm	6,6 cm	354-8222
450 cc	12,7 cm	10,8 cm	7,0 cm	354-8223
550 cc	13,5 cm	11,7 cm	7,4 cm	354-8224
650 cc	14,6 cm	12,6 cm	7,6 cm	354-8225
800 cc	15,6 cm	13,3 cm	8,0 cm	354-8226

CPX™ 4 Breast Tissue Expander, Tall Height, Style 8300*



Volume	Width	Height	Projection	Catalog #
250 cc	10,1 cm	10,7 cm	5,6 cm	354-8321
350 cc	11,3 cm	11,8 cm	6,0 cm	354-8322
450 cc	12,3 cm	12,9 cm	6,5 cm	354-8323
550 cc	13,2 cm	13,8 cm	6,9 cm	354-8324
650 cc	14,0 cm	14,6 cm	7,3 cm	354-8325
750 cc	14,6 cm	15,3 cm	7,6 cm	354-8326
850 cc	15,4 cm	15,9 cm	7,9 cm	354-8327

*Includes CENTERSCOPE™ Magnetic Detection Device

PLEASE NOTE: CONTOUR PROFILE™ Tissue Expanders are not MRI compatible

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

CONTOUR PROFILE™ Tissue Expanders

**CPX™ 4 Breast Tissue Expander with Suture Tabs
Low Height, Style 9100***

Volume	Width	Height	Projection	Catalog #
250 cc	11,4 cm	8,1 cm	6,1 cm	354-9121
350 cc	12,7 cm	9,4 cm	6,5 cm	354-9122
450 cc	14,0 cm	10,2 cm	7,1 cm	354-9123
550 cc	15,0 cm	10,9 cm	7,4 cm	354-9124
650 cc	15,7 cm	11,2 cm	7,9 cm	354-9125
750 cc	16,5 cm	11,9 cm	8,1 cm	354-9126

**CPX™ 4 Breast Tissue Expander with Suture Tabs
Medium Height, Style 9200***

Volume	Width	Height	Projection	Catalog #
275 cc	10,7 cm	9,3 cm	6,2 cm	354-9221
350 cc	11,7 cm	10,0 cm	6,6 cm	354-9222
450 cc	12,7 cm	10,8 cm	7,0 cm	354-9223
550 cc	13,5 cm	11,7 cm	7,4 cm	354-9224
650 cc	14,6 cm	12,6 cm	7,6 cm	354-9225
800 cc	15,6 cm	13,3 cm	8,0 cm	354-9226

**CPX™ 4 Breast Tissue Expander with Suture Tabs
Tall Height, Style 9300***

Volume	Width	Height	Projection	Catalog #
250 cc	10,1 cm	10,7 cm	5,6 cm	354-9321
350 cc	11,3 cm	11,8 cm	6,0 cm	354-9322
450 cc	12,3 cm	12,9 cm	6,5 cm	354-9323
550 cc	13,2 cm	13,8 cm	6,9 cm	354-9324
650 cc	14,0 cm	14,6 cm	7,3 cm	354-9325
750 cc	14,6 cm	15,3 cm	7,6 cm	354-9326
850 cc	15,4 cm	15,9 cm	7,9 cm	354-9327

*Includes CENTERSCOPE™ Magnetic Detection Device

PLEASE NOTE: CONTOUR PROFILE™ Tissue Expanders are not MRI compatible


Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

TISSUE EXPANDERS

Other Tissue Expanders

Smooth Round Tissue Expanders




Device Volume	Diameter	Projection	Catalog #
400 cc	11,3 cm	5,9 cm	350-4305M
550 cc	13,4 cm	5,8 cm	350-4307M
700 cc	14,6 cm	6,7 cm	350-4304M
850 cc	15,0 cm	7,2 cm	350-4308M
1000 cc	15,5 cm	7,8 cm	350-4306M

Smooth Elliptical Tissue Expanders




Device Volume	Length	Width	Nominal Projection	TEMPORARY OVEREXPANSION		Catalog #
				Projection	Max. Volume	
25 cc	6,0 cm	4,0 cm	2,3 cm	4,5 cm	75 cc	350-5301M
50 cc	8,0 cm	5,0 cm	2,5 cm	5,3 cm	150 cc	350-5302M
75 cc	10,0 cm	6,0 cm	2,1 cm	6,0 cm	225 cc	350-5303M
150 cc	12,0 cm	7,0 cm	3,2 cm	7,2 cm	450 cc	350-5304M
250 cc	14,0 cm	8,0 cm	4,5 cm	9,0 cm	750 cc	350-5305M
125 cc	15,0 cm	5,0 cm	3,3 cm	6,5 cm	375 cc	350-5306M

Smooth Crescent Tissue Expanders



Device Volume	Length	Width	Projection	Catalog #
50 cc	9,0 cm	5,7 cm	2,7 cm	350-4312M
100 cc	11,0 cm	6,9 cm	3,7 cm	350-4313M
200 cc	14,6 cm	7,5 cm	4,3 cm	350-4314M
300 cc	16,5 cm	10,3 cm	4,7 cm	350-4315M
400 cc	18,2 cm	11,4 cm	5,2 cm	350-4316M
500 cc	19,6 cm	12,4 cm	5,5 cm	350-4317M

Smooth Rectangle Tissue Expanders



Device Volume	Length	Width	Projection	Catalog #
50 cc	6,9 cm	3,2 cm	3,2 cm	350-4309M
100 cc	6,6 cm	5,1 cm	5,6 cm	350-4301M
250 cc	10,4 cm	7,3 cm	5,6 cm	350-4302M
400 cc	10,6 cm	9,3 cm	6,6 cm	350-4303M
550 cc	13,7 cm	7,1 cm	7,7 cm	350-4310M
700 cc	15,0 cm	8,0 cm	7,9 cm	350-4311M

Sterile Saline Breast Implant Sizers*

Round Saline Breast Implant Sizers, Moderate Profile



NOMINAL FILL			MAXIMUM FILL			Catalog #
Volume	Diameter	Projection	Volume	Diameter	Projection	
125 cc	9,5 cm	3,0 cm	150 cc	9,1 cm	3,4 cm	351-125SZ
175 cc	10,6 cm	3,3 cm	200 cc	10,0 cm	3,9 cm	351-175SZ
225 cc	11,5 cm	3,5 cm	250 cc	10,9 cm	4,0 cm	351-225SZ
275 cc	12,3 cm	3,7 cm	300 cc	11,7 cm	4,2 cm	351-275SZ
325 cc	13,0 cm	3,8 cm	375 cc	12,3 cm	4,5 cm	351-325SZ
375 cc	13,6 cm	4,0 cm	425 cc	12,8 cm	4,8 cm	351-375SZ
425 cc	14,2 cm	4,1 cm	475 cc	13,5 cm	4,9 cm	351-425SZ
475 cc	14,8 cm	4,2 cm	525 cc	14,2 cm	4,9 cm	351-475SZ

Round Saline Breast Implant Sizers, Moderate Plus Profile



NOMINAL FILL			MAXIMUM FILL			Catalog #
Volume	Diameter	Projection	Volume	Diameter	Projection	
175 cc	9,5 cm	3,5 cm	210 cc	9,3 cm	4,2 cm	351-2175SZ
200 cc	10,0 cm	3,7 cm	240 cc	9,7 cm	4,4 cm	351-2200SZ
225 cc	10,4 cm	3,8 cm	270 cc	10,1 cm	4,6 cm	351-2225SZ
250 cc	10,8 cm	4,0 cm	300 cc	10,5 cm	4,7 cm	351-2250SZ
275 cc	11,0 cm	4,1 cm	330 cc	10,8 cm	4,9 cm	351-2275SZ
300 cc	11,5 cm	4,3 cm	360 cc	11,2 cm	5,0 cm	351-2300SZ
325 cc	11,9 cm	4,4 cm	390 cc	11,5 cm	5,2 cm	351-2325SZ
350 cc	12,1 cm	4,5 cm	420 cc	11,7 cm	5,3 cm	351-2350SZ
375 cc	12,3 cm	4,6 cm	450 cc	12,0 cm	5,4 cm	351-2375SZ
400 cc	12,6 cm	4,7 cm	480 cc	12,3 cm	5,5 cm	351-2400SZ
425 cc	12,9 cm	4,8 cm	510 cc	12,5 cm	5,6 cm	351-2425SZ
450 cc	13,0 cm	4,9 cm	540 cc	12,8 cm	5,7 cm	351-2450SZ
475 cc	13,3 cm	5,0 cm	570 cc	13,0 cm	5,8 cm	351-2475SZ
500 cc	13,6 cm	5,1 cm	600 cc	13,2 cm	5,9 cm	351-2500SZ
550 cc	14,0 cm	5,3 cm	660 cc	13,7 cm	6,1 cm	351-2550SZ
600 cc	14,5 cm	5,5 cm	720 cc	14,1 cm	6,3 cm	351-2600SZ
650 cc	14,8 cm	5,6 cm	780 cc	14,4 cm	6,5 cm	351-2650SZ
700 cc	15,2 cm	5,8 cm	840 cc	14,8 cm	6,6 cm	351-2700SZ
750 cc	15,6 cm	5,9 cm	900 cc	15,1 cm	6,8 cm	351-2750SZ
800 cc	16,1 cm	6,0 cm	960 cc	15,5 cm	7,0 cm	351-2800SZ

Note: Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above

SALINE BREAST IMPLANTS

Sterile Saline Implant Sizers[†] CONT'D

Round Saline Breast Implant Sizers, High Profile



NOMINAL FILL			MAXIMUM FILL			Catalog #
Volume	Diameter	Projection	Volume	Diameter	Projection	
190 cc	9,3 cm	4,1 cm	225 cc	9,1 cm	5,1 cm	351-3190SZ
230 cc	10,0 cm	4,3 cm	275 cc	9,8 cm	5,5 cm	351-3230SZ
270 cc	10,4 cm	4,6 cm	325 cc	10,2 cm	5,8 cm	351-3270SZ
290 cc	10,8 cm	4,7 cm	350 cc	10,5 cm	5,9 cm	351-3290SZ
310 cc	11,0 cm	4,8 cm	375 cc	10,7 cm	6,0 cm	351-3310SZ
330 cc	11,3 cm	4,8 cm	400 cc	11,0 cm	6,2 cm	351-3330SZ
380 cc	11,7 cm	5,2 cm	450 cc	11,4 cm	6,4 cm	351-3380SZ
420 cc	12,0 cm	5,4 cm	500 cc	11,7 cm	6,7 cm	351-3420SZ
500 cc	12,8 cm	5,6 cm	600 cc	12,4 cm	7,1 cm	351-3500SZ
630 cc	13,8 cm	5,9 cm	750 cc	13,4 cm	7,4 cm	351-3630SZ

CONTOUR PROFILE™ Saline Breast Implant Sizers, High Profile



Volume	Width	Height	Projection	Catalog #
275 + 25 cc	11,5 cm	9,5 cm	5,1 cm	351-0711SZ
350 + 50 cc	12,3 cm	10,5 cm	5,3 cm	351-0712SZ
450 + 50 cc	13,2 cm	11,0 cm	6,1 cm	351-0713SZ
550 + 50 cc	14,0 cm	11,9 cm	6,4 cm	351-0714SZ
650 + 75 cc	15,0 cm	12,7 cm	6,6 cm	351-0715SZ

Accessories

Injection Port Detector

Description	Catalog #
CENTERSCOPE™ Magnetic Detection Device	350-4402

Injection Dome

Description	Catalog #
Dome Pack Set	350-DOMPK
Can be used with SPECTRUM™ Adjustable Saline Breast Implant, Round BECKER™ Expander/Breast Implant, and Smooth Tissue Expanders.	

Important Safety Information

MENTOR® Breast Implants are indicated for breast augmentation, in women who are at least 18 years old, or for breast reconstruction. Breast implant surgery should not be performed in women with active infection anywhere in their body, with existing cancer or pre-cancer of their breast(s) who have not received adequate treatment for those conditions or who are pregnant or nursing.

There are risks associated with breast implant surgery. Breast implants are not lifetime devices and breast implantation is not necessarily a one-time surgery.

The most common complications with MENTOR® MemoryGel™ Breast Implants include re-operation, implant removal, capsular contracture, asymmetry, and breast pain. A lower risk of complication is implant rupture, which is most often silent. The health consequences of a ruptured silicone gel-filled breast implant have not been fully established. Screenings such as mammography, MRI, or ultrasound are recommended after initial implant surgery to assist in detecting implant rupture.

The most common complications with MENTOR® Saline-Filled Implants include re-operation, implant removal, capsular contracture, wrinkling, deflation, asymmetry, and breast pain.

MENTOR® CPX™4 Breast Tissue Expanders can be utilized for breast reconstruction after mastectomy, correction of an underdeveloped breast, scar revision and tissue defect procedures. These expanders are intended for temporary subcutaneous or submuscular implantation; they should be used within a time frame determined by the physician to achieve the clinically desired degree of tissue expansion.*

CONTOUR PROFILE™ Tissue Expanders are devices that contain magnetic injection domes and are NOT MRI compatible. Do not use the CONTOUR PROFILE™ Tissue Expander in patients where an MRI may be needed. DO NOT use the CONTOUR PROFILE™ Tissue Expander in patients that have a previously implanted device that could be affected by a magnetic field. The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas or where severe surgical reduction of the area has previously been performed.

Your patient needs to be informed and understand the risks and benefits of breast implants, and provided with an opportunity to consult with you prior to deciding on surgery.

For detailed indications, contraindications, warning and precautions associated with the use of all MENTOR® Implantable Devices, please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided at www.mentorwllc.eu.

*** Note:**

CPX™4 Expander is available with two separate Indications for Use. One that is limited to use not more than six months, and the other with the indications shown here. This is dependent upon which CPX™ 4 Expander product code is ordered. Please refer to the Product Insert Data Sheet to determine which product codes are at your customer's availability.

[illegible]

MENTOR WORLDWIDE LLC

Santa Barbara, CA 93111 USA

Tel: +1 805 879 6000

Fax: +1 805 967 7108

www.mentorwwllc.com

Operations in Europe, Middle East & Africa

Italy

Johnson & Johnson Medical S.p.a.

Tel: +39 06 911 94 500

Fax: +39 02 644 10 563

Germany

Mentor Deutschland GmbH

Tel: +49 811 600 50 0

Fax: +49 811 600 50 11

Switzerland

Johnson & Johnson AG -

Franchise Mentor

Tel : +41 800 830 085

Fax: +41 800 820 085

France

ETHICON France, S.A.S.

Tel: +33 1 55 00 22 45

Fax: +33 1 55 00 28 11

Spain & Portugal

Johnson & Johnson, S.A.

Tel: +34 91 562 2700

Fax: +34 91 728 3367

United Kingdom & Ireland

Johnson & Johnson Medical Ltd

Tel: +44 (0)800 973 082

Fax: +44 (0)134 486 4122

Austria

Johnson & Johnson Medical

Phone: +43 1 360 250

Email: jnjmedicalaustria@its.jnj.com

Benelux

Johnson & Johnson Medical BV

Phone: +32 2 481 7400

Email: jnjcsmed@its.jnj.com

Greece

Johnson & Johnson Hellas Commercial
and Industrial S.A.

Tel: +30 210 6875751

Fax: + 30 210 6850308

Israel

Johnson & Johnson Medical Israel

Div. of J-C Healthcare

Tel: +97 2 9959 1111

Sweden, Denmark, Norway,

Finland & Iceland

Johnson & Johnson AB

Tel: +46 8 626 2200

Fax: +46 8 626 2230

Croatia, Macedonia, Serbia, Slovenia, Bulgaria & Romania

Johnson & Johnson S.E. d.o.o.

Tel: +386 1 40 11 800

Fax: +386 1 40 11 801

Czech Republic, Slovak Republic & Hungary

Johnson & Johnson s.r.o.

Tel: +420 227 012 371

Fax: +420 227 012 300

Lithuania, Latvia & Estonia

Johnson & Johnson, UAB

Tel: +370 52 63 61 62

Fax: +370 52 75 50 08

Poland

Johnson & Johnson Poland Sp. z o.o.

Tel: +48 22 237 83 05

Fax: +48 22 237 83 53

Russia & CIS

Johnson & Johnson LLC

Tel: + 7 495 580 7777

Fax: +7 495 580 7878

Turkey

Johnson & Johnson Medical Turkey

Tel: +90 216 538 21 21

Fax: +90 216 538 21 30

Middle East, Egypt & Pakistan

Johnson & Johnson Middle East Inc.

Tel: +971 4 429 72 00

Fax: +971 4 299 43 30

Maghreb & Sub Saharan Africa

Johnson & Johnson Professional Export

Tel: + 44 (0)1506 594 500

Fax: + 44 (0)1506 594 731

South Africa

Johnson & Johnson Medical (Pty.) Ltd.

Tel: +27 11 265 1135

Fax: +27 11 265 1114